





















Competitive Cross-Reference Guide











NiTi Files

| ProTaper®* | | ISO | Channels™ | | | |
|---|---------|-----|-----------|----|---------|--|
|  | purple | S1 | | G1 | purple |  |
| | white | S2 | | G2 | white | |
|  | yellow | F1 | 20 | C1 | yellow |  |
|  | red | F2 | 25 | C2 | red |  |
|  | blue | F3 | 30 | C3 | blue |  |
|  | black | F4 | 40 | C4 | black |  |
|  | yellow | F5 | 50 | C5 | yellow |  |
| | no band | SX | | GX | no band | |

Paper Points/Gutta Percha

| ProTaper®* | | ISO | Channels™ | | | |
|---|--------|-----|-----------|----|--------|---|
|  | yellow | F1 | 20 | C1 | yellow |  |
|  | red | F2 | 25 | C2 | red |  |
|  | blue | F3 | 30 | C3 | blue |  |
|  | black | F4 | 40 | C4 | black |  |
|  | yellow | F5 | 50 | C5 | yellow |  |

Obturators

| ProTaper®* | | ISO | Channels™ | | | |
|---|--------|-----|-----------|----|--------|---|
|  | yellow | F1 | 20 | C1 | yellow |  |
|  | red | F2 | 25 | C2 | red |  |
|  | blue | F3 | 30 | C3 | blue |  |
|  | black | F4 | 40 | C4 | black |  |
|  | yellow | F5 | 50 | C5 | yellow |  |

*ProTaper® is a registered trademark of Dentsply International Inc.